



Fitzgerald

www.fitzgerald-fii.com

Plasma & Blood Proteins



Plasma & Blood Proteins

Plasma or Blood proteins (often termed Serum proteins) are proteins which circulate in blood plasma. Plasma proteins have a variety of different functions including, transport of biological components, regulation of acellular activity, and also play key roles in the immune system. Other plasma proteins constitute enzymes and inhibitors, complement components, and kinin precursors. All plasma proteins are synthesized by the liver except gamma Globulins.

Serum albumin accounts for 55% of blood proteins and plays a major role in maintenance of osmotic pressure of plasma which is key to transporting lipids and steroid hormones. Globulins, which make up about 38% of plasma proteins, are responsible for the transport of ions, hormones, and lipids which function in immune response. Fibrinogen which is an essential component in blood clotting comprises approximately 7% of plasma proteins, while the rest (1%) is made up of regulatory proteins such as enzymes, proenzymes, and hormones.

Current research around blood plasma proteins centers on performing proteomics analyses of serum and plasma in the search for biomarkers. Fitzgerald Industries lists hundreds of antibodies against plasma proteins and purified preparations of plasma proteins many of which are highlighted in this brochure. Please visit our website www.fitzgerald-fii.com for further details.

alpha 1 Antichymotrypsin and alpha 1 Antitrypsin

ANTIBODIES

Catalog No.	Product Name	Host	Clone	Applications
10-A08A	alpha 1 Antichymotrypsin antibody	Mouse	M251002	ELISA
10-A08B	alpha 1 Antichymotrypsin antibody	Mouse	M251001	ELISA
70R-AG003	alpha 1 Antichymotrypsin antibody	Goat	-	User optimized
70R-7584	alpha 1 Antichymotrypsin antibody	Sheep	-	User optimized
10R-1801	alpha 1 Antitrypsin antibody	Mouse	2B12	ELISA, WB
10-A09A	alpha 1 Antitrypsin antibody	Mouse	M2103113	ELISA
20-AG51	alpha 1 Antitrypsin antibody	Goat	-	User optimized
70R-7560	alpha 1 Antitrypsin antibody	Goat	-	User optimized
60R-AG003bt	alpha 1 Antitrypsin antibody (biotin)	Goat	-	ELISA, WB
70R-10678	alpha 1 Antitrypsin antibody	Rabbit	-	ICC, IHC
20R-AR009	alpha 1 Antitrypsin antibody	Rabbit	-	IHC-F, IHC-P
20-AR51	alpha 1 Antitrypsin antibody	Rabbit	-	IHC-F, IHC-P
20C-CR1021SP	alpha 1 Antitrypsin antibody	Sheep	-	ODD, Immunodiffusion, WB

PROTEINS

Catalog No.	Product Name	Protein Type	Grade & Purity	Applications
30-1036	alpha 1 Antichymotrypsin protein	Native	> 90% pure	User optimized
30R-AA002	alpha 1 Antichymotrypsin protein	Native	> 90% pure	User optimized
30-AC47	alpha 1 Antichymotrypsin protein	Native	> 95% pure	User optimized
30-1037	alpha 1 Antitrypsin protein	Native	> 95% pure	User optimized
30R-AA009	alpha 1 Antitrypsin protein	Native	> 96% pure	User optimized
30-AA23	alpha 1 Antitrypsin protein	Native	> 95% pure	User optimized
80R-1217	alpha 1 Antichymotrypsin protein (His tag)	Recombinant	> 95% pure	SDS-PAGE

alpha 1 and alpha 2 Antiplasmin

ANTIBODIES

Catalog No.	Product Name	Host	Clone	Applications
70R-10579	alpha 1 Antiplasmin antibody	Goat	-	User optimized
20R-1362	alpha 1 Antiplasmin antibody	Goat	-	User optimized
70R-10580	alpha 1 Antiplasmin antibody	Sheep	-	User optimized
20R-1356	alpha 1 Antiplasmin antibody	Sheep	-	User optimized
60R-1008	alpha 1 Antiplasmin antibody (HRP)	Sheep	-	User optimized
10-1504	alpha 2 Antiplasmin antibody	Mouse	M59758	ELISA (Detection)
10-1503	alpha 2 Antiplasmin antibody	Mouse	M59757	ELISA (Detection)
10-1502	alpha 2 Antiplasmin antibody	Mouse	M59756	ELISA (Coating)
10R-A123a	alpha 2 Antiplasmin antibody	Mouse	14AP	ELISA
10R-A123b	alpha 2 Antiplasmin antibody	Mouse	3AP	ELISA
10R-A123c	alpha 2 Antiplasmin antibody	Mouse	2AP	ELISA
20R-1393	alpha 2 Antiplasmin antibody	Goat	-	User optimized
60R-1075	alpha 2 Antiplasmin antibody (HRP)	Goat	-	User optimized
70R-10607	alpha 2 Antiplasmin antibody	Goat	-	User optimized
70R-13338	alpha 2 Antiplasmin antibody	Rabbit	-	IF, ICC, IHC-P, WB
20R-1400	alpha 2 Antiplasmin antibody	Sheep	-	User optimized

PROTEINS

Catalog No.	Product Name	Protein Type	Grade & Purity	Applications
30C-CP3004	alpha 2 Antiplasmin protein	Native	> 95% pure	User optimized

alpha 1 Acid Glycoprotein

ANTIBODIES

Catalog No.	Product Name	Host	Clone	Applications
20-OG20	alpha 1 Acid Glycoprotein antibody	Goat	-	User optimized
20C-CR6001G	alpha 1 Acid Glycoprotein antibody	Goat	-	User optimized

PROTEINS

Catalog No.	Product Name	Protein Type	Grade & Purity	Applications
30-1034	alpha 1 Acid Glycoprotein protein	Native	> 96% pure	User optimized
30R-AA003	alpha 1 Acid Glycoprotein protein	Native	> 95% pure	User optimized
30C-CP3001U	alpha 1 Acid Glycoprotein protein	Native	> 95% pure	User optimized
30-AA01	alpha 1 Acid Glycoprotein protein	Native	> 95% pure	User optimized
30-1035	alpha 1 Acid Glycoprotein protein (Mouse)	Native	> 95% pure	User optimized

alpha 1 Microglobulin

ANTIBODIES

Catalog No.	Product Name	Host	Clone	Applications
10R-8508	alpha 1 Microglobulin antibody	Mouse	10A12	ELISA, WB

PROTEINS

Catalog No.	Product Name	Protein Type	Grade & Purity	Applications
30-1039	alpha 1 Microglobulin protein	Native	> 98% pure	User optimized
30-AM06	alpha 1 Microglobulin protein	Native	> 90% pure	User optimized
30R-AM037-10	alpha 1 Microglobulin protein	Native	> 96% pure	User optimized
30-AM05	alpha 1 Microglobulin protein	Native	> 98% pure	User optimized
30C-CP1021	alpha 1 Microglobulin protein	Native	> 95% pure	User optimized

alpha 2 Macroglobulin

ANTIBODIES

Catalog No.	Product Name	Host	Clone	Applications
10R-7499	alpha 2 Macroglobulin antibody	Mouse	3D1	IP, ELISA
10R-7498	alpha 2 Macroglobulin antibody	Mouse	9A3	IP, ELISA
10R-7497	alpha 2 Macroglobulin antibody	Mouse	2D9	WB, ELISA
10C-CR2005M1	alpha 2 Macroglobulin antibody	Mouse	3127485	User optimized
70R-10581	alpha 2 Macroglobulin antibody	Goat	-	User optimized
20R-1357	alpha 2 Macroglobulin antibody	Goat	-	User optimized
20R-MG002	alpha 2 Macroglobulin antibody	Goat	-	ELISA, IHC, WB
20-MG33	alpha 2 Macroglobulin antibody	Goat	-	User optimized
60R-1009	alpha 2 Macroglobulin antibody (HRP)	Goat	-	User optimized
20C-CR2005SP	alpha 2 Macroglobulin antibody	Sheep	-	ODD, Immunodiffusion, WB

PROTEINS

Catalog No.	Product Name	Protein Type	Grade & Purity	Applications
30-1040	alpha 2 Macroglobulin protein	Native	> 96% pure	User optimized
30R-AM036	alpha 2 Macroglobulin protein	Native	> 96% pure	User optimized
30C-CP3005	alpha 2 Macroglobulin protein	Native	> 95% pure	User optimized
30-1041	alpha 2 Macroglobulin protein (Rat)	Native	> 90% pure	User optimized
30-1042	alpha Macroglobulin protein	Native	> 95% pure	User optimized

beta 2 Microglobulin

ANTIBODIES

Catalog No.	Product Name	Host	Clone	Applications
10-1561	beta 2 Microglobulin antibody	Mouse	M64647	ELISA, Latex Immunoturbidimetric Assay
10-1560	beta 2 Microglobulin antibody	Mouse	M64646	ELISA, Latex Immunoturbidimetric Assay
10R-M154c	beta 2 Microglobulin antibody	Mouse	B2	ELISA
10R-M154b	beta 2 Microglobulin antibody	Mouse	C21	User optimized
20-MG30	beta 2 Microglobulin antibody	Goat	-	User optimized
20-MG30IG	beta 2 Microglobulin antibody	Goat	-	User optimized
70R-7564	beta 2 Microglobulin antibody	Goat	-	User optimized
70R-7565	beta 2 Microglobulin antibody (HRP)	Goat	-	User optimized
70R-5904	beta 2 Microglobulin antibody	Rabbit	-	WB
20R-MR013	beta 2 Microglobulin antibody	Rabbit	-	ELISA, IHC, WB
60R-MR005BT	beta 2 Microglobulin antibody (biotin)	Rabbit	-	ELISA, IHC, WB
60R-MR009HRP	beta 2 Microglobulin antibody (HRP)	Rabbit	-	ELISA, IHC, WB
20-MS30	beta 2 Microglobulin antibody	Sheep	-	User optimized
20C-CR2010SP	beta 2 Microglobulin antibody	Sheep	-	ODD, Immunodiffusion, WB

PROTEINS

Catalog No.	Product Name	Protein Type	Grade & Purity	Applications
30-1359	beta 2 Microglobulin protein	Native	> 95% pure	ELISA, Latex Immunoturbidimetric Assay
30R-3114	beta 2 Microglobulin protein	Native	> 95% pure	ELISA, SDS-PAGE, WB
30C-CP1022P	beta 2 Microglobulin protein	Native	> 40% pure	User optimized
30C-CP1022	beta 2 Microglobulin protein	Native	> 95% pure	User optimized
30R-AM001	beta 2 Microglobulin protein	Native	> 95% pure	User optimized
30C-CP1022U	beta 2 Microglobulin protein	Native	> 98% pure	User optimized
30R-3197	beta 2 Microglobulin protein	Recombinant	> 90% pure	ELISA, SDS-PAGE, WB
80R-3423	beta 2 Microglobulin protein (His tag)	Recombinant	> 95% pure	ELISA, SDS-PAGE, WB
80R-1386	beta 2 Microglobulin protein (His tag)	Recombinant	> 95% pure	SDS-PAGE

Ceruloplasmin

ANTIBODIES

Catalog No.	Product Name	Host	Clone	Applications
10R-10409	Ceruloplasmin antibody	Mouse	M78779	ELISA, WB
10R-8442	Ceruloplasmin antibody	Mouse	3B11	ELISA, WB, IP, IHC-P
20-S6006G000-S4	Ceruloplasmin antibody	Goat	-	User optimized
70R-15288	Ceruloplasmin antibody	Rabbit	-	ELISA, IHC, WB
60R-1701	Ceruloplasmin antibody (biotin)	Rabbit	-	ELISA, IHC, WB
60R-1700	Ceruloplasmin antibody (FITC)	Rabbit	-	ELISA, IHC, WB
60R-1699	Ceruloplasmin antibody (HRP)	Rabbit	-	ELISA, IHC, WB
20C-CR6010SP	Ceruloplasmin antibody	Sheep	-	User optimized

PROTEINS

Catalog No.	Product Name	Protein Type	Grade & Purity	Applications
30C-CP1010	Ceruloplasmin protein	Native	> 95% pure	User optimized
30-1075	Ceruloplasmin protein (Mouse)	Native	> 90% pure	User optimized
30-1076	Ceruloplasmin protein (Rat)	Native	> 90% pure	User optimized

C-Reactive Protein (CRP)

ANTIBODIES

Catalog No.	Product Name	Host	Clone	Applications
10R-10264	CRP antibody	Mouse	M78627	WB, IP
10R-8426	CRP antibody	Mouse	26D7	ELISA, WB
10-C189A	CRP antibody	Mouse	M12619CR2	ELISA (Detection), Immunoassays, IHC, WB
10-C189B	CRP antibody	Mouse	M12619CR1	ELISA (Coating), Immunoassays, IHC, WB
10-1004	CRP antibody	Mouse	7A5F21E5	User optimized
10-7894	CRP antibody	Mouse	M18612CRP2	ELISA (Detection), Immunoassays, IHC, WB
10-7893	CRP antibody	Mouse	M18612CRP1	ELISA (Coating), Immunoassays, IHC, WB
10R-C189a	CRP antibody	Mouse	C6	ELISA, IHC, ELISA, TBA
10-C33A	CRP antibody	Mouse	M701289	ELISA (Coating), Immunoassays
10-C33C	CRP antibody	Mouse	M701288	ELISA (Detection), IHC, TBA
10-C33AS	CRP antibody	Mouse	M701289	ELISA (Detection), IHC, TBA
10-C33B	CRP antibody	Mouse	M7111422	ELISA, IHC, ELISA, TBA
10R-C189b	CRP antibody	Mouse	C2	User optimized
10R-3002	CRP antibody	Mouse	KT39	ELISA
10R-1042	CRP antibody	Mouse	5A9	ELISA, WB
10C-CR2015M6	CRP antibody	Mouse	3127455	User optimized
10C-CR2015M5	CRP antibody	Mouse	2127455	User optimized
70R-15116	CRP antibody	Goat	-	ELISA, IHC, WB
20-1002	CRP antibody	Goat	-	User optimized
70-B9007GA00-A0	CRP antibody	Goat	-	User optimized
20-S6011G000-S4	CRP antibody	Goat	-	User optimized
20-S3911G001-V0	CRP antibody	Goat	-	User optimized
20-S3911G000-V0	CRP antibody	Goat	-	User optimized
20-S3911G000-S4	CRP antibody	Goat	-	User optimized
20-S3903GND9-D0	CRP antibody	Goat	-	User optimized
20-B9007GD00-D0	CRP antibody	Goat	-	User optimized
20-B9007G000-W0	CRP antibody	Goat	-	User optimized
20-S3903GND1-D0	CRP antibody	Goat	-	User optimized
70-XG20	CRP antibody	Goat	-	User optimized
70-XG21	CRP antibody	Goat	-	User optimized
20-CG15	CRP antibody	Goat	-	User optimized
20C-CR2015GAP	CRP antibody	Goat	-	User optimized
70-XG20S	CRP antibody	Goat	-	User optimized
70R-13710	CRP antibody	Rabbit	-	WB
70R-4527	CRP antibody	Rabbit	-	WB
20-1003	CRP antibody	Rabbit	-	User optimized
70R-10663	CRP antibody	Rabbit	-	IHC, ICC, WB
20-1262	CRP antibody	Sheep	-	ELISA
20C-CR2015SP	CRP antibody	Sheep	-	ODD, Immunodiffusion, WB
20C-CR2015S	CRP antibody	Sheep	-	User optimized
60R-1194	CRP antibody (biotin)	Goat	-	ELISA, IHC, WB
60R-1195	CRP antibody (FITC)	Goat	-	ELISA, IHC, WB
60R-1196	CRP antibody (HRP)	Goat	-	ELISA, IHC, WB

PROTEINS

Catalog No.	Product Name	Protein Type	Grade & Purity	Applications
30-1094	CRP protein	Native	> 98% pure	
30-AC10	CRP protein	Native	Partially pure	Control
30-S3911H000-T0	CRP protein	Native	-	User optimized
30R-AC067	CRP protein	Native	> 99% pure	User optimized
30C-CP1000U	CRP protein	Native	Highly pure	Control
30R-AC001x	CRP protein	Native	Highly pure	Control
30-AC05PP	CRP protein	Native	> 80% pure	
30-AC05S	CRP protein	Native	Highly pure	
30-AC05AF	CRP protein (Azide free)	Native	> 95% pure	User optimized
30-1093	CRP protein (Canine)	Native	> 96% pure	
30-1095	CRP protein (Rat)	Native	> 96% pure	
30R-2081	CRP protein	Recombinant	> 90% pure	SDS-PAGE
30-1092	CRP protein	Recombinant	> 95% pure	
30-AC07	CRP protein	Recombinant	-	User optimized
30R-2740	CRP protein	Recombinant	> 95% pure	
30R-2389	CRP protein	Recombinant	> 95% pure	
30R-2984	CRP protein	Recombinant	> 95% pure	

Fibrinogen

ANTIBODIES

Catalog No.	Product Name	Host	Clone	Applications
10R-8419	Fibrinogen antibody	Mouse	5C5	ELISA, WB, IHC-P
10-F107C	Fibrinogen antibody	Mouse	M12619FB2	ELISA (Detection), Immunoassays, WB
10-F107B	Fibrinogen antibody	Mouse	M12619FB1	ELISA (Coating/Detection), Immunoassays, WB
10-7903	Fibrinogen antibody	Mouse	M18612FB	ELISA (Coating), Immunoassays, WB
10R-3003	Fibrinogen antibody	Rat	KT9	ELISA
10R-F107C	Fibrinogen antibody	Mouse	27C8	ELISA, WB
10R-F107B	Fibrinogen antibody	Mouse	1F3	ELISA, WB
70R-7578	Fibrinogen antibody	Goat	-	User optimized
20C-CR6029RP	Fibrinogen antibody	Rabbit	-	ELISA, RID, Electrophoresis
70R-10582	Fibrinogen antibody	Sheep	-	User optimized
70R-10575	Fibrinogen antibody	Sheep	-	User optimized
20R-1366	Fibrinogen antibody	Sheep	-	User optimized
20R-1354	Fibrinogen antibody	Sheep	-	User optimized
20C-CR6029SP	Fibrinogen antibody	Sheep	-	ODD, Immunodiffusion, WB
61R-1023	Fibrinogen antibody (HRP)	Rat	KT10	ELISA
60R-1013	Fibrinogen antibody (HRP)	Sheep	-	User optimized
60R-1012	Fibrinogen antibody (HRP)	Sheep	-	User optimized
60R-1011	Fibrinogen antibody (HRP)	Sheep	-	User optimized
60R-FS003hrp	Fibrinogen antibody (HRP)	Sheep	-	User optimized
20R-1420	Fibrinogen antibody (IgG fraction)	Goat	-	User optimized
70R-5273	Fibrinogen alpha antibody	Rabbit	-	WB
70R-1553	Fibrinogen alpha antibody	Rabbit	-	WB, IHC

PROTEINS

Catalog No.	Product Name	Protein Type	Grade & Purity	Applications
30C-CP6131P	Fibrinogen protein	Native	> 40% pure	User optimized
30R-AF013	Fibrinogen protein	Native	> 95% pure	User optimized
30-AF27	Fibrinogen protein	Native	> 95% pure	User optimized
30C-CP6131U	Fibrinogen protein	Native	> 98% pure	User optimized
30R-3109	Fibrinogen protein (Bovine)	Native	> 95% pure	ELISA, SDS-PAGE, WB

Globulin

ANTIBODIES

Catalog No.	Product Name	Host	Clone	Applications
20-S4111G000-V0	Gc Globulin antibody	Goat	-	User optimized
20-S4111G000-S2	Globulin antibody	Goat	-	User optimized
20-B9010G000-S0	Gamma globulin antibody	Goat	-	User optimized

PROTEINS

Catalog No.	Product Name	Protein Type	Grade & Purity	Applications
30-1111	Gamma Globulin protein	Native	> 99% pure	User optimized
30-AB44	Bovine gamma Globulin	Native	> 99% pure	User optimized
30R-GS001	Pig gamma Globulin	Native	-	User optimized
30-AI36	Human gamma Globulin protein	Native	-	User optimized

Haptoglobin

ANTIBODIES

Catalog No.	Product Name	Host	Clone	Applications
10R-10318	Haptoglobin antibody	Mouse	M78688	ELISA, WB
10R-8440	Haptoglobin antibody	Mouse	26E11	ELISA, WB
10C-CR2114M1	Haptoglobin antibody	Mouse	637400	User optimized
20-1012	Haptoglobin antibody	Goat	-	Nephelometric Assay, Latex Immunoturbidimetric Assay
70R-14136	Haptoglobin antibody	Rabbit	-	WB, IHC-P
70R-13514	Haptoglobin antibody	Rabbit	-	IF, ICC, IHC-P, WB
70R-12755	Haptoglobin antibody	Rabbit	-	IF, ICC, IHC-P, WB
20-S7401G000-S4	Haptoglobin antibody	Goat	-	User optimized
20-S1011G000-V0	Haptoglobin antibody	Goat	-	User optimized
70R-5459	Haptoglobin antibody	Rabbit	-	WB
70R-5419	Haptoglobin antibody	Rabbit	-	WB
70C-CR2114RAP	Haptoglobin antibody	Rabbit	-	ELISA, IHC, WB
20-HG10	Haptoglobin antibody	Goat	-	User optimized
20C-CR2114G	Haptoglobin antibody	Goat	-	User optimized
20C-CR2114GP	Haptoglobin antibody	Goat	-	User optimized
70C-CR1338RAP	Haptoglobin antibody	Rabbit	-	User optimized
70C-CR1337GAP	Haptoglobin antibody	Goat	-	User optimized
70C-CR2114GAP	Haptoglobin antibody	Goat	-	User optimized
20C-CR2114SP	Haptoglobin antibody	Sheep	-	ODD, Immunodiffusion, WB
70R-7558	Haptoglobin antibody	Goat	-	User optimized

PROTEINS

Catalog No.	Product Name	Protein Type	Grade & Purity	Applications
30-1126	Haptoglobin protein	Native	> 98% pure	User optimized
30C-CP3018	Haptoglobin protein	Native	> 95% pure	User optimized
30-AH02	Haptoglobin protein	Native	> 95% pure	User optimized
30-1124	Haptoglobin protein (Bovine)	Native	> 90% pure	User optimized
30-1125	Haptoglobin protein (Canine)	Native	> 95% pure	User optimized
30-1127	Haptoglobin protein (Mouse)	Native	> 90% pure	User optimized
30-1128	Haptoglobin protein (Porcine)	Native	> 90% pure	User optimized
30-1129	Haptoglobin protein (Rat)	Native	> 96% pure	User optimized

Plasminogen

ANTIBODIES

Catalog No.	Product Name	Host	Clone	Applications
10-7945	Plasminogen antibody	Mouse	M53470	ELISA (Coating), Immunoassays
10C-CR6065M1	Plasminogen antibody	Mouse	357500	ELISA (Detection), IHC
10R-2516	Plasminogen antibody	Mouse	10 VPA	ELISA
10R-2515	Plasminogen antibody	Mouse	10 VPA	ELISA
10R-2514	Plasminogen antibody	Mouse	9 VPA	ELISA
10R-P132a	Plasminogen antibody	Mouse	4D2	ELISA, IHC-F, IHC-P
10R-8473	Plasminogen antibody	Mouse	44A1	ELISA, IP
10R-8472	Plasminogen antibody	Mouse	9A1	ELISA, IP, IHC-P
10R-8453	Plasminogen antibody	Mouse	10A1	ELISA, WB, IHC-P
20-S4611G000-S4	Plasminogen antibody	Goat	-	User optimized
20-S4611G000-V0	Plasminogen antibody	Goat	-	User optimized
20-S4611G000-S0	Plasminogen antibody	Goat	-	User optimized
20-S4603GND1-D0	Plasminogen antibody	Goat	-	User optimized
70R-10577	Plasminogen antibody	Goat	-	User optimized
20R-1359	Plasminogen antibody	Goat	-	User optimized
20R-PG001	Plasminogen antibody	Goat	-	User optimized
70R-7562	Plasminogen antibody	Goat	-	User optimized
70R-5384	Plasminogen antibody	Rabbit	-	WB
20-PR29	Plasminogen antibody	Rabbit	-	User optimized
20C-CR8063SP	Plasminogen antibody	Sheep	-	ELISA, RIA, WB
70R-10578	Plasminogen antibody	Sheep	-	User optimized
20R-1360	Plasminogen antibody	Sheep	-	User optimized
60R-1007	Plasminogen antibody (HRP)	Goat	-	User optimized

PROTEINS

Catalog No.	Product Name	Protein Type	Grade & Purity	Applications
30-AP27	Plasminogen protein	Native	> 95% pure	Human

Prothrombin

ANTIBODIES

Catalog No.	Product Name	Host	Clone	Applications
10R-8454	Prothrombin precursor antibody	Mouse	28H5	ELISA, WB
70R-10583	Prothrombin factor II antibody	Sheep	-	User optimized
20R-1358	Prothrombin factor II antibody	Sheep	-	User optimized
60R-1044	Prothrombin factor II antibody (HRP)	Sheep	-	User optimized
70R-10602	Prothrombin fragment 1 antibody	Sheep	-	User optimized
70R-10587	Prothrombin fragment 2 antibody	Sheep	-	User optimized

PROTEINS

Catalog No.	Product Name	Protein Type	Grade & Purity	Applications
30C-CP3048U	Prothrombin protein	Native	> 98% pure	User optimized
30C-CP3049U	Prothrombin protein (Bovine)	Native	> 95% pure	User optimized
30-1204	Prothrombin protein (Porcine)	Native	> 95% pure	User optimized

Transferrin

ANTIBODIES

Catalog No.	Product Name	Host	Clone	Applications
10R-T118a	Transferrin antibody	Mouse	2A2	ELISA, IF
10-T34D	Transferrin antibody	Mouse	M10021522	ELISA (Coating)
10-T34E	Transferrin antibody	Mouse	M10021523	ELISA (Detection)
10R-10575	Transferrin antibody	Mouse	M79068	WB
10-T118A	Transferrin antibody	Mouse	M12619TF1	ELISA (Detection), Immunoassays
10-T34C	Transferrin antibody	Mouse	M10021521	ELISA (Coating)
10-T34A	Transferrin antibody	Mouse	M1120701	User optimized
10R-8524	Transferrin antibody	Mouse	14C10	ELISA, WB, IP
10R-8523	Transferrin antibody	Mouse	3C11	WB
20-1020	Transferrin antibody	Goat	-	Nephelometric Assay, Latex Immunoturbidimetric Assay
70C-CR6072GAP	Transferrin antibody	Goat	-	ELISA, ICC, IHC, WB
70R-7559	Transferrin antibody	Goat	-	User optimized
20-S6007G000-S4	Transferrin antibody	Goat	-	User optimized
20-S5493G000-S4	Transferrin antibody	Goat	-	User optimized
20-S1611G000-V0	Transferrin antibody	Goat	-	User optimized
20-S1611G000-S4	Transferrin antibody	Goat	-	User optimized
20-S1603GND1-D0	Transferrin antibody	Goat	-	User optimized
20-D3913GN01-D5	Transferrin antibody	Goat	-	User optimized
20-B9014G000-S0	Transferrin antibody	Goat	-	User optimized
20-TG43	Transferrin antibody	Goat	-	User optimized
20C-CR6072G	Transferrin antibody	Goat	-	User optimized
20-TG43S	Transferrin antibody	Goat	-	User optimized
70R-15486	Transferrin antibody	Rabbit	-	ELISA, IHC, WB
70R-15371	Transferrin antibody	Rabbit	-	ELISA, IHC, WB
20R-2963	Transferrin antibody	Rabbit	-	WB, IP
70R-10681	Transferrin antibody	Rabbit	-	IP, ELISA, IHC
20C-CR6072RP	Transferrin antibody	Rabbit	-	ELISA, WB
20-TR43	Transferrin antibody	Rabbit	-	User optimized
20R-TR038	Transferrin antibody	Rabbit	-	User optimized
20R-TR047	Transferrin antibody	Rabbit	-	User optimized
20R-TR046	Transferrin antibody	Rabbit	-	User optimized
20R-TR045	Transferrin antibody	Rabbit	-	User optimized
60R-2277	Transferrin antibody (biotin)	Rabbit	-	ELISA, IHC, WB
60R-1932	Transferrin antibody (biotin)	Rabbit	-	ELISA, IHC, WB
60R-TR002bt	Transferrin antibody (biotin)	Rabbit	-	ELISA, IHC, WB
60R-2276	Transferrin antibody (FITC)	Rabbit	-	ELISA, IHC, WB
60R-1931	Transferrin antibody (FITC)	Rabbit	-	ELISA, IHC, WB
60R-TR009FT	Transferrin antibody (FITC)	Rabbit	-	User optimized
60R-2275	Transferrin antibody (HRP)	Rabbit	-	ELISA, IHC, WB
60R-1930	Transferrin antibody (HRP)	Rabbit	-	ELISA, IHC, WB

PROTEINS

Catalog No.	Product Name	Protein Type	Grade & Purity	Applications
30-1222	Transferrin Apo protein	Native	> 99% pure	User optimized
30-1223	Transferrin Apo protein (Rat)	Native	> 98% pure	User optimized
30-1224	Transferrin Holo protein	Native	> 96% pure	User optimized
30R-3085	Transferrin protein	Native	> 90% pure	ELISA, SDS-PAGE, WB
30-AT31	Transferrin protein	Native	> 95% pure	User optimized
30C-CP4060	Transferrin protein	Native	> 95% pure	User optimized
30-1225	Transferrin protein	Recombinant	> 98% pure	User optimized
30R-AT019	Transferrin protein (Apo)	Native	> 98% pure	User optimized
30R-3089	Transferrin protein (Bovine)	Native	> 90% pure	ELISA, SDS-PAGE, WB
30R-AT020	Transferrin protein (Holo)	Native	> 98% pure	User optimized
30R-AT021	Transferrin protein (Iron-rich)	Native	> 99% pure	User optimized
31C-CH1026	Transferrin protein (Mouse)	Native	-	ELISA, WB, IHC
30R-3092	Transferrin protein (Rabbit)	Native	> 90% pure	ELISA, SDS-PAGE, WB
30R-3100	Transferrin protein (Rat)	Native	> 90% pure	ELISA, SDS-PAGE, WB
30C-CH1216	Transferrin protein (Rat)	Native	-	ELISA, WB, IHC
30R-3111	Transferrin protein (Sheep)	Native	> 95% pure	ELISA, SDS-PAGE, WB



Fitzgerald

Website www.fitzgerald-fii.com
E-mail antibodies@fitzgerald-fii.com
Toll Free (USA) 800.370.2222
Tel 978.371.6446
Fax 978.371.2266



scan code for more info