




Fitzgerald

Alzheimer's Disease



Click on any category below to view the relevant products

**Monoclonal
Antibodies**

**Polyclonal
Antibodies**


Proteins

**Normal Brain
Lysates**

Research Kits

Alzheimer's Disease: Monoclonal Antibodies

Catalog No.	Product Name	Host	Clone	Size
10R-2281	160 kDa Neurofilament antibody	Mouse	NF403	1 ml
10R-3050	160/210 kDa Neurofilament antibody	Mouse	N2184	500 ul
10R-N119a	200 kDa Neurofilament antibody	Mouse	RT97	100 ug
10R-2280	200 kDa Neurofilament antibody	Mouse	NF402	1 ml
10R-2279	70/200 kDa Neurofilament antibody	Mouse	2F11	1 ml
10R-8054	ABPP antibody	Mouse	C4D10	100 ug
10R-8055	ABPP antibody	Mouse	A10E2	100 ug
10R-1019	alpha Synuclein antibody	Mouse	5C2	100 ul
10R-1837	alpha Synuclein antibody	Mouse	2B2D1, 2B2A11	100 ul
10R-2292	Amyloid A antibody	Mouse	Reu-86.2	1 ml
10R-A126b	Amyloid P antibody	Mouse	5.4A	100 ug
10R-A126a	Amyloid P antibody	Mouse	5.4D.3B	100 ug
10R-A136a	ApoE antibody	Mouse	E3	50 ug
10R-2370	ApoE antibody	Mouse	1B5-E1	50 ug
10R-3015	ApoE antibody	Mouse	5B5	200 ug
10R-A127a	APP antibody	Mouse	40.10	1 ml
10R-1097	APP antibody	Mouse	J4H9	100 ul
10R-2016	APP antibody	Mouse	1B11F3	100 ul
10R-3016	APP antibody	Mouse	10D1	100 ug
10R-7541	beta Amyloid antibody	Mouse	AP-5	100 ug
10R-B113a	beta Amyloid antibody	Mouse	6F/3D	1 ml
10R-A124a	beta Amyloid antibody	Mouse	22C11	50 ug
10R-3011	beta Amyloid antibody	Mouse	12B2	50 ug
10R-3012	beta Amyloid antibody	Mouse	1A10	50 ug
10R-3013	beta Amyloid antibody	Mouse	82E1	50 ug
10R-7917	beta Amyloid antibody	Mouse	AP-5	100 ug
10R-3014	beta Amyloid antibody (biotin)	Mouse	82E1	50 ug
10R-2702	beta Tubulin 3 antibody	Mouse	2E9	100 ul



Immunofluorescent staining of beta Tubulin 3 antibody (10R-2702)


Immunofluorescence analysis of PANC-1 cells using beta tubulin 3 antibody (green). Blue: DRAQ5 fluorescent DNA dye. Red: Actin filaments have been labeled with Alexa Fluor-555 phalloidin.

Alzheimer's Disease: Monoclonal Antibodies

Catalog No.	Product Name	Host	Clone	Size
10R-D117a	Drebrin antibody	Mouse	M2F6	100 ug
10R-D117b	Drebrin antibody	Mouse	MX823	250 ul
10R-2460	Drebrin antibody	Mouse	MX823	5 ml
10R-G105a	GFAP antibody	Mouse	GF2.24	50 ug
10R-2792	GFAP antibody	Mouse	6A6	100 ul
10R-2910	GFAP antibody	Mouse	5C8	100 ul
10R-3036	GFAP antibody	Mouse	G2584	1 ml
10R-7290	GFAP antibody	Mouse	GA-8	100 ug
10R-7859	GFAP antibody	Mouse	GA-8	100 ug
10R-T102b	phospho Tau antibody	Mouse	AT8	100 ug
10R-T102b	phospho Tau antibody	Mouse	AT8	50 ug
10R-3064	Presenilin 1 antibody	Mouse	17C2	200 ug
10R-7419	PTEN antibody	Mouse	30G12	100 ul
10R-3017	sAPP alpha antibody	Mouse	2B3	100 ug
10R-2852	SORL1 antibody	Mouse	3B6B11	100 ul
10R-8037	Tau 13 antibody	Mouse	B11E8	100 ul
10R-T101a	Tau antibody	Mouse	TAU-5	200 ug
10R-T101a	Tau antibody	Mouse	TAU-5	100 ug
10R-3069	Tau antibody	Mouse	9A1	100 ug
10R-3070	Tau antibody	Mouse	1A1	100 ug
10R-3071	Tau antibody	Mouse	C5	100 ug
10R-3072	Tau antibody	Mouse	2B11	100 ug
10-B9091MD01-MD	Tau antibody	Mouse	B9091-Tau	1 mg
10R-7337	Tau antibody	Mouse	TAU-93	100 ug
10R-7903	Tau antibody	Mouse	TAU-93	100 ug

Alzheimer's Disease: Polyclonal Antibodies

Catalog No.	Product Name	Host	Size
70R-10632	160 kDa Neurofilament antibody	Chicken	200 ul
70R-10630	200 kDa Neurofilament antibody	Chicken	200 ul
70R-10631	68 kDa Neurofilament antibody	Chicken	1 ml
20R-1060	AKT antibody (Phospho Ser473)	Rabbit	100 ul
20R-1058	AKT1 antibody	Rabbit	100 ul
70R-10493	AKT1 antibody	Rabbit	50 ug
33R-8182	AKT1 Blocking Peptide	>70R-10493<	100 ug
70R-3068	AKT1S1 antibody	Rabbit	50 ug
33R-4654	AKT1S1 Blocking Peptide	>70R-3068<	100 ug
70R-2220	AKTIP antibody	Rabbit	50 ug
33R-6833	AKTIP Blocking Peptide	>70R-2220<	100 ug
70R-SG007	alpha Synuclein antibody	Goat	100 ug
20R-2019	alpha Synuclein antibody	Rabbit	50 ug
20R-2091	alpha Synuclein antibody	Rabbit	50 ug
20R-2113	alpha Synuclein antibody	Rabbit	50 ug
20R-2114	alpha Synuclein antibody	Rabbit	50 ug
20R-2348	alpha Synuclein antibody	Rabbit	50 ug
20R-2406	alpha Synuclein antibody	Rabbit	50 ug
20R-2432	alpha Synuclein antibody	Rabbit	50 ug
20R-2433	alpha Synuclein antibody	Rabbit	50 ug
20R-2824	alpha Synuclein antibody	Rabbit	50 ug
20R-2825	alpha Synuclein antibody	Rabbit	50 ug
70R-7799	alpha Synuclein antibody	Rabbit	100 ug
20-SG70	alpha Synuclein antibody	Goat	100 ul
20R-SG003	alpha Synuclein antibody	Goat	250 ul
20R-2821	alpha Synuclein antibody (Phospho-Ser129)	Rabbit	50 ug
20R-2822	alpha Synuclein antibody (Phospho-Tyr133)	Rabbit	50 ug
20R-2823	alpha Synuclein antibody (Phospho-Tyr136)	Rabbit	50 ug
20C-CR1325G	Amyloid P antibody	Goat	1 ml
70R-10555	APBA1 antibody	Rabbit	50 ug
33R-7273	APBA1 Blocking Peptide	>70R-10555<	100 ug
70R-4577	APBB3 antibody	Rabbit	50 ug
33R-8950	APBB3 Blocking Peptide	>70R-4577<	100 ug
70C-CR9016GAP	ApoE antibody	Goat	100 ug
20C-CR9016G	ApoE antibody	Goat	100 ul




Immunohistochemical staining using APBA1 antibody (70R-10555)

APBA1 antibody was used for immunohistochemistry at a concentration of 4-8 ug/ml to stain Epithelial cells of intestinal villus (arrows) in Human Intestine. Magnification is at 400X.

Alzheimer's Disease: Polyclonal Antibodies

Catalog No.	Product Name	Host	Size
20-AG45	ApoE antibody	Goat	1 ml
70R-5264	ApoE antibody	Rabbit	50 ug
33R-4716	ApoE Blocking Peptide	>70R-5264<	100 ug
70R-7639	ApoE antibody	Rabbit	100 ug
70R-AR011	ApoE3 antibody	Rabbit	50 ug
20-AG53	ApoJ antibody	Goat	500 ug
70R-12314	ApoSAA antibody	Rabbit	100 ug
20C-CR1088GP	APP antibody	Goat	1 mg
20R-AG004	APP antibody	Goat	250 ul
20-AG70	APP antibody	Goat	100 ul
20-AG71	APP antibody	Goat	100 ul
20R-AG005	APP antibody	Goat	250 ul
20-AG72	APP antibody	Goat	100 ul
70R-6194	APP antibody	Rabbit	50 ug
33R-8067	APP Blocking Peptide	>70R-6194<	100 ug
70R-6195	APP antibody	Rabbit	50 ug
33R-2119	APP Blocking Peptide	>70R-6195<	100 ug
70R-6196	APP antibody	Rabbit	50 ug
33R-2264	APP Blocking Peptide	>70R-6196<	100 ug
20R-1034	APP antibody	Rabbit	100 ul
70R-7640	APP antibody	Rabbit	100 ug
70R-7641	APP antibody	Rabbit	100 ug
70R-7642	APP antibody	Rabbit	100 ug
70R-7643	APP antibody	Rabbit	100 ug
70R-9609	APP antibody	Rabbit	50 ug
33R-5529	APP Blocking Peptide	>70R-9609<	100 ug
70R-11129	APP antibody	Rabbit	50 ug
70R-12464	APP antibody	Rabbit	100 ug
70R-13905	APP antibody	Rabbit	100 ug
70R-11128	APP antibody (Phospho-668)	Rabbit	50 ug
70R-2199	APPBP2 antibody	Rabbit	50 ug
33R-7479	APPBP2 Blocking Peptide	>70R-2199<	100 ug




Immunohistochemical staining using ApoE antibody (70R-5264)

ApoE antibody was used for immunohistochemistry at a concentration of 4-8 ug/ml. Magnification is at 400X.

Alzheimer's Disease: Polyclonal Antibodies

Catalog No.	Product Name	Host	Size
70R-BG005	BACE antibody	Goat	100 ug
20R-BR014	BACE1 antibody	Rabbit	50 ug
20R-BR015	BACE1 antibody	Rabbit	50 ug
70R-7276	BACE1 antibody	Rabbit	50 ug
33R-3496	BACE1 Blocking Peptide	>70R-7276<	100 ug
70R-7645	BACE1 antibody	Rabbit	100 ug
70R-7646	BACE1 antibody	Rabbit	100 ug
70R-8662	BACE1 antibody	Rabbit	50 ug
70R-6559	BACE2 antibody	Rabbit	50 ug
33R-8612	BACE2 Blocking Peptide	>70R-6559<	100 ug
70R-6662	BACE2 antibody	Rabbit	50 ug
33R-6958	BACE2 Blocking Peptide	>70R-6662<	100 ug
70R-6209	BDNF antibody	Rabbit	50 ug
33R-2813	BDNF Blocking Peptide	>70R-6209<	100 ug
70R-10695	beta Amyloid A4 antibody	Rabbit	50 ug
70R-10694	beta Amyloid A4 antibody (Phospho-Thr743/668)	Rabbit	50 ug
70R-BR019	beta Amyloid antibody	Rabbit	100 ug
70R-AR020	beta Amyloid antibody	Rabbit	100 ug
20R-AR011	beta Amyloid antibody	Rabbit	250 ul
70R-7633	beta Amyloid antibody	Rabbit	50 ug
70R-7634	beta Amyloid antibody	Rabbit	50 ug
70R-7635	beta Amyloid antibody	Rabbit	50 ug
70R-7636	beta Amyloid antibody	Rabbit	100 ug
70R-7637	beta Amyloid antibody	Rabbit	50 ug
70R-10624	beta Amyloid antibody	Chicken	1 ml
20R-2839	beta Amyloid antibody	Rabbit	50 ul
20R-1000	beta Synuclein antibody	Rabbit	100 ul
70R-10639	beta Tubulin 3 antibody	Chicken	1 ml



Immunohistochemical staining using BACE1 antibody (70R-7276)


BACE1 antibody was used for immunohistochemistry at a concentration of 4-8 ug/ml. Magnification is at 400X.

Alzheimer's Disease: Polyclonal Antibodies

Catalog No.	Product Name	Host	Size
70R-10534	CHAT antibody	Rabbit	50 ug
33R-8800	CHAT Blocking Peptide	>70R-10534<	100 ug
70R-10625	CHAT antibody	Chicken	1 ml
70R-CR039	CLACP antibody	Rabbit	100 ug
70R-CR040	CLACP antibody	Rabbit	100 ug
70R-CR041	CLACP antibody	Rabbit	100 ug
70R-7682	Drebrin A antibody	Rabbit	50 ug
20R-DP004	Drebrin antibody	Guinea Pig	100 ul
20R-DP003	Drebrin antibody	Guinea Pig	100 ul
70R-GR031	GFAP antibody	Rabbit	100 ug
70R-GR030	GFAP antibody	Rabbit	100 ug
20C-CR8013R	GFAP antibody	Rabbit	100 ul
20R-GR005	GFAP antibody	Rabbit	250 ul
20R-GR004	GFAP antibody	Rabbit	1 ml
20R-GR004	GFAP antibody	Rabbit	500 ul
70R-10640	GFAP antibody	Chicken	200 ul
70R-10675	GFAP antibody	Rabbit	500 ul
70R-10675	GFAP antibody	Rabbit	1 ml
20R-2732	GFAP antibody	Rabbit	50 ug
70R-10705	GFAP antibody	Rabbit	100 ug
70R-13832	GFAP antibody	Rabbit	100 ug
75R-1041	GFAP antibody (Prediluted for IHC)	Rabbit	7 ml
70R-6020	MAPT antibody	Rabbit	50 ug
33R-6638	MAPT Blocking Peptide	>70R-6020<	100 ug
70R-6021	MAPT antibody	Rabbit	50 ug
33R-7914	MAPT Blocking Peptide	>70R-6021<	100 ug
70R-14021	MAPT antibody	Rabbit	100 ug
20R-BR016	Mouse Brain antibody	Rabbit	2 ml
20R-BR019	Mouse Brain antibody	Rabbit	2 ml
20R-BR017	Mouse Brain antibody	Rabbit	1 ml
20R-BR018	Mouse Brain antibody	Rabbit	1 ml
60R-BR001FT	Mouse Brain antibody (FITC)	Rabbit	2 mg
60R-BR002FT	Mouse Brain antibody (FITC)	Rabbit	2 mg
70R-MG008	MSR Type 1 antibody	Goat	100 ug
70R-NR016	Neprilysin antibody	Rabbit	100 ug
20R-NR019	NFT antibody	Rabbit	100 ug
20-NR65	NFT antibody	Rabbit	250 ul
20R-NR004	NFT antibody	Rabbit	250 ul
70R-2589	PADI2 antibody	Rabbit	50 ug
33R-7920	PADI2 Blocking Peptide	>70R-2589<	100 ug
70R-12472	Phospho-NR2B antibody	Rabbit	100 ul
70R-12468	Phospho-Tau antibody	Rabbit	100 ul

Alzheimer's Disease: Polyclonal Antibodies

Catalog No.	Product Name	Host	Size
20-PG71	Presenilin 1 antibody	Goat	100 ul
20C-CR1220GP	Presenilin 1 antibody	Goat	1 mg
70R-6264	Presenilin 1 antibody	Rabbit	50 ug
33R-9260	Presenilin 1 Blocking Peptide	>70R-6264<	100 ug
70R-PG003	Presenilin 2 antibody	Goat	100 ug
20-PG72	Presenilin 2 antibody	Goat	100 ul
70R-5685	Presenilin 2 antibody	Rabbit	50 ug
33R-9902	Presenilin 2 Blocking Peptide	>70R-5685<	100 ug
70R-13901	Presenilin 2 antibody	Rabbit	100 ug
70R-9773	PSENEN antibody	Rabbit	50 ug
33R-8714	PSENEN Blocking Peptide	>70R-9773<	100 ug
20-PS57	PTEN antibody	Sheep	250 ug
20R-1922	PTEN antibody	Rabbit	50 ug
20R-2237	PTEN antibody	Rabbit	50 ug
20R-2238	PTEN antibody	Rabbit	50 ug
70R-10976	PTEN antibody	Rabbit	50 ug
70R-10977	PTEN antibody	Rabbit	50 ug
70R-10978	PTEN antibody	Rabbit	50 ug
70R-10975	PTEN antibody (Phospho-Ser370)	Rabbit	50 ug
70R-10973	PTEN antibody (Phospho-Ser380)	Rabbit	50 ug
70R-10974	PTEN antibody (Phospho-Ser380/Thr382/Thr383)	Rabbit	50 ug
70R-RR002	Rattin antibody	Rabbit	100 ug
33R-6638	MAPT Blocking Peptide	>70R-6020<	100 ug
70R-6021	MAPT antibody	Rabbit	50 ug
70R-6965	RTN1 antibody	Rabbit	50 ug
33R-2384	RTN1 Blocking Peptide	>70R-6965<	100 ug
70R-1758	RTN2 antibody	Rabbit	100 ug



Immunohistochemical staining using Presenilin 2 antibody (70R-5685)

Presenilin 2 antibody was used for immunohistochemistry at a concentration of 4-8 ug/ml to stain Epithelial cells of renal tubule (arrows) in Human Kidney. Magnification is at 400X.

Alzheimer's Disease: Polyclonal Antibodies

Catalog No.	Product Name	Host	Size
33R-6064	RTN2 Blocking Peptide	>70R-1758<	100 ug
70R-7137	RTN4 antibody	Rabbit	50 ug
33R-7833	RTN4 Blocking Peptide	>70R-7137<	100 ug
70R-7139	RTN4 antibody	Rabbit	50 ug
33R-3051	RTN4 Blocking Peptide	>70R-7139<	100 ug
20-SG20	SAA antibody	Goat	1 ml
70R-7644	sAPP beta Wild Type antibody	Rabbit	100 ug
70R-SG006	Synphilin 1 antibody	Goat	100 ug
20R-SG002	Synphilin 1 antibody	Goat	250 ul
70R-TG010	Tau antibody	Goat	100 ug
20R-TG001	Tau antibody	Goat	250 ul
70R-7800	Tau antibody	Rabbit	100 ug
70R-10637	Tau antibody	Chicken	1 ml
20R-1958	Tau antibody	Rabbit	50 ug
20R-1959	Tau antibody	Rabbit	50 ug
20R-1963	Tau antibody	Rabbit	50 ug
20R-1964	Tau antibody	Rabbit	50 ug
20R-1965	Tau antibody	Rabbit	50 ug
20R-1966	Tau antibody	Rabbit	50 ug
20R-1967	Tau antibody	Rabbit	50 ug
20R-1968	Tau antibody	Rabbit	50 ug
20R-1969	Tau antibody	Rabbit	50 ug
20R-2092	Tau antibody	Rabbit	50 ug
20R-2102	Tau antibody	Rabbit	50 ug
20R-2182	Tau antibody	Rabbit	50 ug
20R-2272	Tau antibody	Rabbit	50 ug
20R-2273	Tau antibody	Rabbit	50 ug
20R-2275	Tau antibody	Rabbit	50 ug
20R-2276	Tau antibody	Rabbit	50 ug
20R-2277	Tau antibody	Rabbit	50 ug
20R-2278	Tau antibody	Rabbit	50 ug
20R-2279	Tau antibody	Rabbit	50 ug
20R-2280	Tau antibody	Rabbit	50 ug
20R-2281	Tau antibody	Rabbit	50 ug
20R-2407	Tau antibody	Rabbit	50 ug
20R-2420	Tau antibody	Rabbit	50 ug
20R-2692	Tau antibody	Rabbit	50 ug
70R-11904	Tau antibody	Rabbit	100 ug
33R-10799	Tau Blocking Peptide	>70R-11904<	50 ug

Alzheimer's Disease: Proteins

Catalog No.	Product Name	Species	Source	Size
30R-2112	160 kDa Neurofilament protein	Bovine	Bovine spinal cord	250 ug
30R-2132	160 kDa Neurofilament protein	Bovine	Bovine spinal cord	100 ug
30R-AN002	160 kDa Neurofilament protein	Bovine	Bovine spinal cord	250 ul
30R-2113	200 kDa Neurofilament protein	Bovine	Bovine spinal cord	250 ug
30R-2133	200 kDa Neurofilament protein	Bovine	Bovine spinal cord	100 ug
30R-AN037	200 kDa Neurofilament protein	Bovine	Bovine spinal cord	250 ug
30R-2111	68 kDa Neurofilament protein	Bovine	Bovine spinal cord	250 ug
30R-2131	68 kDa Neurofilament protein	Bovine	Bovine spinal cord	100 ug
30R-AN003	68 kDa Neurofilament protein	Bovine	Bovine spinal cord	250 ug
30R-1009	alpha Synuclein 160 protein	Human	Recombinant	100 ug
30R-1013	alpha Synuclein 195 protein	Human	Recombinant	100 ug
30R-1007	alpha Synuclein 61140 protein	Human	Recombinant	100 ug
30R-1008	alpha Synuclein 96140 protein	Human	Recombinant	100 ug
30R-1002	alpha Synuclein A30P protein	Human	Recombinant	100 ug
30R-1004	alpha Synuclein A30P/A53T protein	Human	Recombinant	100 ug
30R-1003	alpha Synuclein A53T protein	Human	Recombinant	100 ug
30R-1005	alpha Synuclein delta NAC protein	Human	Recombinant	100 ug
30R-1012	alpha Synuclein E46K protein	Human	Recombinant	100 ug
30R-1000	alpha Synuclein protein	Human	Recombinant	100 ug
30R-1001	alpha Synuclein protein	Human	Recombinant	100 ug
30C-CP1104	Amyloid P protein	Human	Human plasma	250 ug
30-AA30	ApoE protein	Human	Human plasma	100 ug
30R-AA019	ApoE2 protein	Human	Recombinant	500 ug
30R-2386	ApoE2 protein	Human	Recombinant	100 ug
30R-2386	ApoE2 protein	Human	Recombinant	50 ug
30R-2387	ApoE2 protein	Human	Recombinant	500 ug
30R-2387	ApoE2 protein	Human	Recombinant	5 mg
30R-2387	ApoE2 protein	Human	Recombinant	1 mg
30R-AA016	ApoE3 protein	Human	Recombinant	500 ug
30R-2381	ApoE3 protein	Human	Recombinant	5 mg
30R-2381	ApoE3 protein	Human	Recombinant	1 mg
30R-AA020	ApoE4 protein	Human	Recombinant	500 ug
30R-AA046	ApoE4 protein	Human	Recombinant	100 ug
30R-2382	ApoE4 protein	Human	Recombinant	5 mg
30R-2382	ApoE4 protein	Human	Recombinant	1 mg
30R-AA017	ApoSAA1 protein	Human	Recombinant	50 ug
30R-AA045	ApoSAA1 protein	Human	Recombinant	50 ug
30R-2420	BDNF protein	Human	Recombinant	10 ug

Alzheimer's Disease: Proteins

Catalog No.	Product Name	Species	Source	Size
30R-2815	beta Secretase protein	Human	Recombinant	20 ug
30R-1010	beta Synuclein protein	Human	Recombinant	100 ug
30R-2110	GFAP protein	Bovine	Bovine spinal cord	250 ug
30R-AG009	GFAP protein	Bovine	Bovine spinal cord	250 ug
80R-1550	KCNIP3 protein (His tag)	Human	Recombinant	100 ug
80R-1468	QPRT protein (His tag)	Human	Recombinant	100 ug
30-1211	Serum Amyloid P protein	Human	Human Serum	250 ug
80R-1408	Tau protein (His tag)	Human	Recombinant	50 ug
80R-1556	Tau protein (His tag)	Human	Recombinant	50 ug
80R-1575	Tau protein (His tag)	Human	Recombinant	50 ug
80R-1787	Tau protein (His tag)	Human	Recombinant	10 ug

Alzheimer's Disease: Normal Human Brain Lysates

Catalog No.	Product Name	Size
30R-AM039	Human Brain Medulla Tissue Lysate	150 ug
30R-AB017	Human Brain Tissue Lysate	150 ug
30R-AC058	Human Cerebellum Tissue Lysate	150 ug
30R-AF	Human Frontal Cortex Tissue Lysate	150 ug
30R-AH	Human Hippocampus Tissue Lysate	150 ug
30R-AM011	Human Midbrain Tissue Lysate	150 ug
30R-AO004	Human Olfactory Lobe Tissue Lysate	150 ug
30R-AP033	Human Pons Tissue Lysate	150 ug
30R-AP034	Human Posterior Cortex Lysate	150 ug
30R-AS030	Human Striatum Tissue Lysate	150 ug
30R-AT053	Human Thalamus Tissue Lysate	150 ug

Alzheimer's Disease: Research Kits

Catalog No.	Product Name
55R-1383	beta Secretase Activity Assay Kit
55R-27418	Human Amyloid beta ELISA kit
55R-27716	Human Amyloid beta ELISA kit
55R-27725	Human Amyloid beta Oligomers ELISA kit
55R-1228	Human APP ELISA Kit
55R-1726	Human APP ELISA Kit
55R-27729	Human beta Amyloid ELISA kit
55R-27717	Human beta Amyloid ELISA kit
55R-27713	Human beta Amyloid ELISA kit
55R-27718	Human beta Amyloid ELISA kit
55R-27714	Human beta Amyloid ELISA kit
55R-27711	Human beta Amyloid ELISA kit
55R-27712	Human beta Amyloid ELISA kit
55R-1198	Human Cystatin C ELISA Kit
55R-RD192072200R	Human GFAP ELISA kit
55R-27769	Human LRG ELISA kit
55R-27734	Human sAPP alpha ELISA kit
55R-27732	Human sAPP betaw ELISA kit
55R-27762	Human/Mouse/Rat alpha 2, 6 Sialyltransferase ELISA kit
55R-27752	Human/Mouse/Rat BACE1 ELISA kit
55R-27731	Human sAPP ELISA kit
55R-1199	Mouse Cystatin C ELISA Kit
55R-1734	Mouse Cystatin C ELISA Kit
55R-27720	Mouse/Rat Amyloid beta ELISA kit
55R-27415	Mouse/Rat sAPP alpha ELISA kit
55R-27408	Rat alpha 2,6 Sialyltransferase ELISA kit
55R-27765	Rat NERC ELISA kit